



OPERATING UNIT







Technical characteristics

- Operating voltage: 24 V
- Air: 5 bar, 50 l/min

- Power supply: 230 v/115 v, 250 W
- ➤ Weight: 115 Kg

CHEST OF DRAWERS' CONFIGURATION



WITHOUT HANDLES

Drawers mounted on "push" slides that allow the opening of the drawer with a slight pressure on the front panel.



WHITE HANDLES


Drawers mounted on soft closing slides and equipped with a newly designed aluminum handle painted in white.





COLORED HANDLES

Drawers mounted on soft closing slides and equipped with a new design aluminum handle painted in a color of your choice from Stern Podia selection.


TYPE OF DRAWERS (T1/T2/T3) AND DOORS (T5/T8):







Possibility to configure the podiatric unit with an additional chest of drawers.





DRAWERS' INTERNAL INSERTS

